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Directions for use of improved Schiötz Tonometer

Before use, the tonometer has to be placed upon the convex test block. The pointer should be on zero. If the pointer has another position, the instrument requires readjustment.

The weight marked 5.5 always remains on the plunger so that the other weights must be added supplementarily to the weight 5.5.

For measuring the tension, the patient is laid in the recumbent position and the cornea anesthetized with two or three drops of 2% Holocain or 1/2% Pontocain. While the patient is instructed to look straight upwards; the lids are retracted gently with one hand, without exerting pressure on the eyeball. The other hand seizes the handle of the tonometer, placing it in a vertical position on the center of the cornea. The handle of the tonometer is lowered to a position midway between the top and foot plate of the cylinder; thereby the instrument will act independently by its own weight. Now the position of the pointer can be read. The tension in mm/Hg may be determined by the table, or by the diagram.

The average normal tension is 15 to 20 mm Hg. A tension of more than 20 mm Hg may already be too high while a tension above 25 mm Hg is considered suspicious and is usually indicative of a pathological condition.

After use, the plunger should be removed from the cylinder. This is done by unscrewing the weight 5.5 g. The plunger, the groove and the foot plate must be cleaned with ether and rubbed dry in order to avoid corrosion.